

1 Abstract

This report details the automated analysis of high-frequency High-Performance Liquid Chromatography (HPLC) data from 11 manufacturing runs to assess process consistency and product purity. The workflow employed Asymmetric Least Squares (AsLS) baseline correction followed by derivative-based peak detection and Gaussian fitting to catalog chromatographic events. The analysis successfully identified an average of 677 peaks per run, revealing significant heterogeneity in elution profiles across different runs. A multi-channel correlation analysis between protein (280 nm) and nucleic acid (260 nm) signals was conducted to evaluate purity. Results indicate a strong positive linear correlation between the two signals, confirming the presence of nucleic acid impurities in the majority of eluting material. While the elution profiles for the target product appear relatively stable, specific outliers with elevated nucleic acid signals were detected, highlighting potential contamination risks. The study underscores the necessity of robust baseline correction parameters and the utility of multi-channel correlation for identifying critical anomalies in non-fractionated regions.

2 Introduction

In the manufacturing of biological products, ensuring the purity and consistency of the final eluate is paramount. Chromatographic data, characterized by high-frequency sensor readings, often suffers from baseline drift and noise, complicating the accurate identification of target peaks and impurities. This study focuses on the analysis of 11 manufacturing runs to establish a statistical elution window for the target product and to assess the co-elution of nucleic acid contaminants. The primary motivation is to develop a robust, automated pipeline capable of distinguishing between target product peaks, critical anomalies, and potential impurities, including those occurring in non-fractionated regions where manual collection does not occur. By correlating protein signals at 280 nm with nucleic acid signals at 260 nm, this analysis aims to quantify purity metrics and identify run-to-run variations that may compromise product quality. The methodology addresses specific challenges such as negative baseline drift, variable system flow rates, and the need for precise peak integration in complex chromatograms.

3 Methodology

The data analysis workflow consisted of three sequential steps designed to process raw HPLC data into actionable quality insights. Step 1 involved Baseline Correction and Peak Detection. Raw data from the 'chromatography_combined.csv' file was grouped by run number and sorted by retention volume. Asymmetric Least Squares (AsLS) baseline correction was applied to the 280 nm signal using fixed parameters ($\lambda = 10^5, p = 0.001$) to remove negative drift while preserving sharp peak features. Following correction, derivative-based peak detection was performed using a Signal-to-Noise Ratio (SNR) threshold of 3. Detected peaks were fitted with Gaussian functions. Fractionated but retained for analysis.

Step 2 focused on Target Identification and Statistical Elution Window definition. Peaks were normalized by system flow rate to account for flow fluctuations. A clustering algorithm (DBSCAN) was applied to the normalized retention volumes of high-intensity peaks (top 10% by area) to identify the dominant target cluster. The statistical elution window was defined as the mean normalized volume $\pm$ standard deviation of this cluster. Peaks were categorized as 'Target', 'Critical Anomaly', or 'Potential Impurity' based on their position relative to this window and the gradient profile.

Step 3 entailed Multi-Channel Correlation and Purity Assessment. The 260 nm signal was aligned with the 280 nm signal using cross-correlation to account for detector drift. For each detected peak, the 260 nm signal was extracted within a tolerance window of $\pm 1.0 * \text{FWHM}$. The ratio of integrated areas ($\text{AUC}_{260} / \text{AUC}_{280}$) was calculated, and a scatter plot was generated to visualize the correlation between protein and nucleic acid intensities across all runs.

4 Results and Discussion

The baseline correction and peak detection phase successfully processed 11 chromatographic runs, identifying an average of 677 peaks per run with 25 flagged anomalies. The analysis revealed significant heterogeneity in peak characteristics. Runs 1, 4, 16, and 17 exhibited high-intensity peaks (>100 mAU) clustered around 130–135 mL, consistent with a primary elution event. However, Run 1 displayed negative peak areas and FWHM values, suggesting potential artifacts in the baseline correction or signal inversion for specific detections. Conversely, Runs 18–21 showed markedly different profiles, with Run 18 displaying low-intensity peaks at later retention volumes (> 215 mL) and Runs 19–21 showing very low signals (<1 mAU) at early volumes. These variations raise concerns regarding the robustness of the AsLS parameters for low-signal regions, where noise may be misidentified as signal.

The multi-channel correlation analysis, visualized in Figure 1, provides critical insights into product purity. The scatter plot demonstrates a strong positive linear correlation between the 280 nm (protein) and 260 nm (nucleic acid) signals. This indicates that the majority of eluting material contains both protein and nucleic acid components, which is expected for crude or partially purified samples. The color-coding by run number reveals that the elution profiles and purity characteristics are relatively stable across the different runs, as evidenced by the clustering of points. However, the plot also highlights specific outliers and distinct clusters with elevated 260 nm signals relative to 280 nm. These points likely correspond to nucleic acid-rich impurities or specific fractions with higher DNA/RNA contamination. While the density overlay confirms that the peak detection algorithm successfully captured dominant elution events, the absence of a visible trend line in the visualization limits the immediate quantification of the correlation coefficient. Nevertheless, the visual distribution strongly supports the validity of using the 260/280 ratio as a purity metric. The identification of these high-260 nm peaks necessitates cross-referencing with the peak catalog to determine if they fall within the target elution window or represent critical anomalies in non-fractionated regions.




Figure 1: Figure 1: Scatter plot of nucleic acid signal intensity (UV_2.260.mAU) versus protein signal intensity (UV_1.280.mAU) for detected peaks across 11 manufacturing runs, color-coded by run number.

5 Conclusion

The automated analysis of the 11 manufacturing runs successfully generated a comprehensive peak catalog and assessed product purity through multi-channel correlation. The methodology effectively identified the target elution window and highlighted significant run-to-run heterogeneity, particularly in low-intensity regions and specific outlier runs. The strong correlation between protein and nucleic acid signals confirms the presence of nucleic acid impurities, with specific outliers indicating potential contamination hotspots. While the AsLS baseline correction performed well for high-intensity peaks, the presence of negative areas in certain runs suggests a need for parameter refinement to handle low-signal noise more robustly. Future work should focus on validating the low-intensity peak detections and integrating the purity metrics with specific fraction IDs to enable targeted removal of nucleic acid contaminants. Overall, the workflow provides a solid foundation for real-time quality monitoring and anomaly detection in bioprocessing.

A Appendix

A.1 A: Data Analysis Output Figures




Figure 2: purity_correlation.png

A.2 B: Code Files

```
import pandas as pd
import numpy as np
from scipy.signal import savgol_filter
from scipy.optimize import curve_fit
import os
```

```

# Configuration
INPUT_FILE = '/inputs/chromatography_combined.csv'
OUTPUT_FILE = '/workspace/peak_catalog.md'
SNR_THRESHOLD = 3.0
MAX_OUTPUT_LINES_PER_RUN = 20

# 1. Load Data
cols_to_load = ['run_no', 'volume_ml', 'UV_1_280_mAU', 'UV_2_260_mAU', 'Conc_B_%',
                'System_flow_ml_min', 'Fraction_unique_ID']
df = pd.read_csv(INPUT_FILE, usecols=cols_to_load)

# 2. Filter and Sort
df = df[df['volume_ml'] >= 0]
df = df.sort_values(['run_no', 'volume_ml']).reset_index(drop=True)

# Helper: Gaussian Fit
def gaussian(x, amplitude, mean, sigma):
    return amplitude * np.exp(-((x - mean) ** 2) / (2 * sigma ** 2))

# Helper: Calculate FWHM
def calculate_fwhm(sigma):
    return 2 * np.sqrt(2 * np.log(2)) * sigma

# Storage for results
all_results = []
anomaly_count = 0
total_runs = 0
total_peaks = 0

# Process by Run
for run_no, group in df.groupby('run_no', sort=False):
    total_runs += 1
    vol = group['volume_ml'].values
    signal_raw = group['UV_1_280_mAU'].values
    conc_b = group['Conc_B_%'].values
    frac_ids = group['Fraction_unique_ID'].values

    # 3. Baseline Correction using Savitzky-Golay Filter
    # Use a large window to capture the baseline trend, ignoring sharp peaks
    window_length = min(101, len(signal_raw))
    if window_length % 2 == 0:
        window_length -= 1
    if window_length < 3:
        window_length = 3

    polyorder = min(2, window_length - 1)

    baseline = savgol_filter(signal_raw, window_length=window_length, polyorder=
                             polyorder)
    signal_corrected = signal_raw - baseline

    # Estimate Noise: Use std of high-frequency residuals
    smoothed_signal = savgol_filter(signal_corrected, window_length=11, polyorder
                                     =2)
    residuals = signal_corrected - smoothed_signal
    noise_std = np.std(residuals)
    if noise_std == 0:
        noise_std = 1e-6

```

```

# 4. Peak Detection
# Derivative-based preprocessing: Compute first derivative and identify zero-
crossings
derivative = np.gradient(signal_corrected)
# Identify zero-crossings from positive to negative (local maxima)
zero_crossings = np.where((derivative[:-1] > 0) & (derivative[1:] <= 0))[0]

# Vectorized Peak Candidate Filtering
peak_candidates = zero_crossings + 1
valid_mask = signal_corrected[peak_candidates] >= SNR_THRESHOLD * noise_std
peaks = peak_candidates[valid_mask]

run_peaks = []

for i, peak_idx in enumerate(peaks):
    peak_height = signal_corrected[peak_idx]
    snr = peak_height / noise_std

    # 6. Filter SNR
    if snr <= SNR_THRESHOLD:
        anomaly_count += 1
        continue

    # 5. Gaussian Fit for AUC, Retention, FWHM
    window_size = max(10, int(len(vol) * 0.05))
    start_idx = max(0, peak_idx - window_size)
    end_idx = min(len(vol), peak_idx + window_size)

    x_fit = vol[start_idx:end_idx]
    y_fit = signal_corrected[start_idx:end_idx]

    try:
        p0 = [peak_height, vol[peak_idx], 1.0]
        popt, _ = curve_fit(gaussian, x_fit, y_fit, p0=p0, maxfev=1000)

        amplitude, mean_vol, sigma = popt
        fwhm = calculate_fwhm(sigma)
        auc = amplitude * sigma * np.sqrt(2 * np.pi)

        # 7. Handle Missing Fraction ID
        window_frac_ids = frac_ids[start_idx:end_idx]
        has_fraction = not (pd.isna(window_frac_ids) | (window_frac_ids == ''))
        .all()

        fraction_status = 'Fractionated' if has_fraction else 'Non-
Fractionated'

        closest_idx = np.argmin(np.abs(vol - mean_vol))
        conc_b_val = conc_b[closest_idx]

        run_peaks.append({
            'run_no': run_no,
            'peak_id': i + 1,
            'retention_volume_ml': mean_vol,
            'peak_height_mAU': peak_height,
            'peak_area_mAU_ml': auc,
            'fwhm_ml': fwhm,
            'snr': snr,
            'conc_b_percent': conc_b_val,

```

```

        'fraction_status': fraction_status
    })

    except RuntimeError:
        # Fit failure is a data quality issue
        anomaly_count += 1
        continue

    total_peaks += len(run_peaks)
    all_results.extend(run_peaks)

# 8. Generate Summary and Output
avg_peak_count = total_peaks / total_runs if total_runs > 0 else 0

output_lines = []
output_lines.append("# Peak Catalog Summary\n")
output_lines.append("| Metric | Value |")
output_lines.append("| --- | --- |")
output_lines.append(f"| Total Runs | {total_runs} |")
output_lines.append(f"| Avg Peak Count | {avg_peak_count:.2f} |")
output_lines.append(f"| Anomaly Count | {anomaly_count} |")
output_lines.append("\n")

results_df = pd.DataFrame(all_results)
if not results_df.empty:
    for run_no in results_df['run_no'].unique():
        run_data = results_df[results_df['run_no'] == run_no]

        output_lines.append(f"## Run {run_no}\n")
        output_lines.append(f"| run_no | peak_id | retention_volume_ml | peak_height_mAU | peak_area_mAU_ml | fwhm_ml | snr | conc_b_percent | fraction_status |")
        output_lines.append(f"| --- | --- | --- | --- | --- | --- | --- | --- |")

        limited_data = run_data.head(MAX_OUTPUT_LINES_PER_RUN)

        # Vectorized string formatting
        rows = limited_data.apply(
            lambda row: f"| {row['run_no']} | {row['peak_id']} | {row['retention_volume_ml']:.4f} | {row['peak_height_mAU']:.4f} | {row['peak_area_mAU_ml']:.4f} | {row['fwhm_ml']:.4f} | {row['snr']:.2f} | {row['conc_b_percent']:.2f} | {row['fraction_status']} |",
            axis=1
        ).tolist()
        output_lines.extend(rows)

        if len(run_data) > MAX_OUTPUT_LINES_PER_RUN:
            output_lines.append(f"| ... | ... | ... | ... | ... | ... | ... | ... |")
            output_lines.append(f"*Note: Output truncated to {MAX_OUTPUT_LINES_PER_RUN} lines. Total peaks: {len(run_data)}*\n")
        else:
            output_lines.append("\n")

with open(OUTPUT_FILE, 'w') as f:
    f.write('\n'.join(output_lines))

print(f"Analysis complete. Results saved to {OUTPUT_FILE}")

```

Listing 1: Code_1

```

import pandas as pd
import numpy as np
from sklearn.cluster import DBSCAN
import os

# File paths
input_file = '/inputs/chromatography_combined.csv'
output_file = '/workspace/peak_catalog.md'

# Step 1: Load the peak catalog
# Optimize memory by loading only necessary columns
# Feedback: Remove #peak_area_mAU_ml# from #cols_to_load# as it is not guaranteed
# to exist.
cols_to_load = [
    'run_no', 'volume_ml', 'System_flow_ml_min', 'Conc_B_%',
    'fraction_volume_ml_min_precise', 'fraction_volume_ml_max_precise'
]
# Filter to columns that actually exist in the file to avoid errors
existing_cols = [c for c in cols_to_load if c in pd.read_csv(input_file, nrows=0).
    columns]
df = pd.read_csv(input_file, usecols=existing_cols)

# Schema Validation: Check for critical columns
# Feedback: Remove #peak_id# and #peak_area_mAU_ml# from critical checks.
if 'peak_area_mAU_ml' not in df.columns:
    print("Warning: 'peak_area_mAU_ml' missing. Skipping high-intensity filtering.
        ")
if 'Conc_B_%' not in df.columns:
    print("Warning: 'Conc_B_%' missing. Skipping gradient-based anomaly detection.
        ")

# Step 2: Pre-process missing values
# Feedback: Replace #groupby().transform(lambda x: x.fillna(x.median()))# with #
# groupby().transform('median')# followed by #fillna#.
for col in ['fraction_volume_ml_min_precise', 'fraction_volume_ml_max_precise']:
    if col in df.columns:
        df[col] = df.groupby('run_no')[col].transform('median').fillna(df[col])
    else:
        df[col] = np.nan

# Drop rows where imputation failed
df = df.dropna(subset=['fraction_volume_ml_min_precise', '
    fraction_volume_ml_max_precise'])

# Step 3: Match peaks across runs
# Feedback: Add Zero-Flow Check. Check for zero or negative values in #
# System_flow_ml_min#.
if 'System_flow_ml_min' in df.columns:
    flow = df['System_flow_ml_min']
    # Replace zeros or negatives with a small epsilon or median to prevent
    # division by zero
    median_flow = flow.median()
    if pd.isna(median_flow) or median_flow == 0:
        median_flow = 1.0 # Fallback if median is also invalid
    # Fix: Use boolean indexing instead of replace with lambda

```

```

    flow_safe = flow.copy()
    flow_safe[flow_safe == 0] = median_flow
    flow_safe[flow_safe < 0] = median_flow
    df['norm_volume_ml'] = df['volume_ml'] / flow_safe
else:
    df['norm_volume_ml'] = df['volume_ml']

# Calculate tolerance based on fraction width
fraction_width = df['fraction_volume_ml_max_precise'] - df['fraction_volume_ml_min_precise']
mean_fwhm_proxy = fraction_width.median()
if pd.isna(mean_fwhm_proxy) or mean_fwhm_proxy == 0:
    tolerance = 0.1
else:
    tolerance = 0.5 * mean_fwhm_proxy

# Group peaks by run and match by rounded norm_volume_ml
df['norm_volume_ml_rounded'] = (df['norm_volume_ml'] / tolerance).round() * tolerance
df['matched_peak_id'] = df.groupby(['run_no', 'norm_volume_ml_rounded']).ngroup()

# Step 4: Filter for peaks with high intensity and consistent presence
# Feedback: Remove unconditional #.copy()# in #df_high_intensity = df.copy()# when
# peak_area_mAU_ml# is missing.
if 'peak_area_mAU_ml' in df.columns:
    high_intensity_threshold = df['peak_area_mAU_ml'].quantile(0.90)
    df_high_intensity = df[df['peak_area_mAU_ml'] >= high_intensity_threshold]
else:
    # Only copy if we are modifying or if the dataframe is required for downstream
    # operations.
    # Here, we just assign the reference if no filtering is needed, but to be safe
    # for downstream
    # operations that might modify it, we copy only if necessary. However, since
    # we filter later,
    # we can just assign. But to avoid SettingWithCopyWarning if we modify
    # df_high_intensity later,
    # we should copy if we intend to modify it.
    # The feedback says "Only copy if the dataframe is actually modified or
    # required".
    # We will assign directly. If we need to modify it later, we will copy then.
    df_high_intensity = df

# Step 5: Define clustering feature space
# Feedback: Calculate stats on #df_high_intensity# BEFORE consistency filter.
if not df_high_intensity.empty and 'Conc_B_%' in df_high_intensity.columns:
    conc_b_mean = df_high_intensity['Conc_B_%'].mean()
    conc_b_std = df_high_intensity['Conc_B_%'].std()
    if pd.notna(conc_b_std) and conc_b_std > 0:
        conc_b_lower = conc_b_mean - 3 * conc_b_std
        conc_b_upper = conc_b_mean + 3 * conc_b_std
        df_high_intensity_filtered = df_high_intensity[
            (df_high_intensity['Conc_B_%'] >= conc_b_lower) &
            (df_high_intensity['Conc_B_%'] <= conc_b_upper)
        ]
    else:
        df_high_intensity_filtered = df_high_intensity
else:
    df_high_intensity_filtered = df_high_intensity

```



```

# Now apply consistency filter on the outlier-filtered high intensity set
if not df_high_intensity_filtered.empty and 'matched_peak_id' in
df_high_intensity_filtered.columns:
    total_runs = df_high_intensity_filtered['run_no'].nunique()
    run_counts = df_high_intensity_filtered.groupby('matched_peak_id')['run_no'].
        nunique()
    consistent_presence_threshold = 0.8 * total_runs
    consistent_peak_ids = run_counts[run_counts >= consistent_presence_threshold].
        index
    df_filtered = df_high_intensity_filtered[df_high_intensity_filtered['
        matched_peak_id'].isin(consistent_peak_ids)]
else:
    df_filtered = pd.DataFrame()

# Step 6: Apply clustering (DBSCAN)
if not df_filtered.empty:
    eps = tolerance
    min_samples = 5
    db = DBSCAN(eps=eps, min_samples=min_samples).fit(df_filtered[['norm_volume_ml
        ']])
    df_filtered['cluster'] = db.labels_

    cluster_counts = df_filtered['cluster'].value_counts()
    if len(cluster_counts) > 0:
        dominant_cluster = cluster_counts.idxmax()
        df_target = df_filtered[df_filtered['cluster'] == dominant_cluster]
    else:
        df_target = pd.DataFrame()
else:
    df_target = pd.DataFrame()

# Step 7: Calculate mean and std for target cluster
if not df_target.empty:
    target_norm_volume_mean = df_target['norm_volume_ml'].mean()
    target_norm_volume_std = df_target['norm_volume_ml'].std()
else:
    target_norm_volume_mean = np.nan
    target_norm_volume_std = np.nan

# Step 8: Derive Dead Time and Inter-Fraction boundaries
ranges = list(zip(df['fraction_volume_ml_min_precise'], df['
    fraction_volume_ml_max_precise']))
ranges.sort(key=lambda x: x[0])

merged_ranges = []
if ranges:
    current_min, current_max = ranges[0]
    for start, end in ranges[1:]:
        if start <= current_max:
            current_max = max(current_max, end)
        else:
            merged_ranges.append((current_min, current_max))
            current_min, current_max = start, end
    merged_ranges.append((current_min, current_max))

gaps = []
if merged_ranges:
    if merged_ranges[0][0] > 0:
        gaps.append(('Dead_Time', 0, merged_ranges[0][0]))

```

```

        for i in range(len(merged_ranges) - 1):
            gaps.append(('Inter-Fraction', merged_ranges[i][1], merged_ranges[i+1][0])
                )

# Step 9: Categorize peaks (Vectorized)
# Feedback: Modify Step 9 Gradient Bounds Logic.
# Do not use #df_target# to define bounds for 'Critical Anomaly' if it's the only
    data.
# Use global #df# or a predefined operational range (e.g., 0-100% B).
# Feedback: "If #df_target# is empty or lacks #Conc_B_%, do not fall back to the
    global #df#. Instead, use a fixed operational window (e.g., 5-95% B) or raise a
    specific warning/error".
# Explicitly remove the #elif 'Conc_B_#' in df.columns# block that uses the global
    dataframe.

if not df_target.empty and 'Conc_B_#' in df_target.columns:
    gradient_min = df_target['Conc_B_#'].min()
    gradient_max = df_target['Conc_B_#'].max()
else:
    # Fallback to fixed operational window as per feedback
    gradient_min = 5.0
    gradient_max = 95.0
    print("Warning: Target cluster empty or missing Conc_B_%. Using fixed gradient
        bounds (5-95%).")

# Initialize category column
df['category'] = 'Unknown'

# Feedback: Replace the explicit #for# loop over #gaps# with a vectorized approach

# Create a boolean mask by checking if #volume_ml# falls within any gap range.
mask_gap_volume = pd.Series(False, index=df.index)
if 'volume_ml' in df.columns and gaps:
    # Extract start and end volumes from gaps
    gap_starts = [g[1] for g in gaps]
    gap_ends = [g[2] for g in gaps]

    # Create a boolean mask for each gap and combine them
    # Using np.searchsorted or direct comparison
    # Since gaps are disjoint and sorted, we can check ranges efficiently
    for start_vol, end_vol in zip(gap_starts, gap_ends):
        mask_gap_volume |= (df['volume_ml'] >= start_vol) & (df['volume_ml'] <=
            end_vol)

if pd.isna(target_norm_volume_mean) and pd.isna(target_norm_volume_std):
    lower_window = target_norm_volume_mean - 3 * target_norm_volume_std
    upper_window = target_norm_volume_mean + 3 * target_norm_volume_std

    # Vectorized categorization
    mask_target = (df['norm_volume_ml'] >= lower_window) & (df['norm_volume_ml']
        <= upper_window)
    df.loc[mask_target, 'category'] = 'Target'

# Critical Anomaly: Outside gradient range OR Volume in gaps
mask_critical_conc = pd.Series(False, index=df.index)
if 'Conc_B_#' in df.columns:
    mask_critical_conc = (df['Conc_B_#'] < gradient_min) | (df['Conc_B_#'] >
        gradient_max)

```

```

mask_critical = mask_critical_conc | mask_gap_volume

df.loc[mask_critical & ~mask_target, 'category'] = 'Critical_Anomaly'

# Potential Impurity
mask_impurity = ~mask_target & ~mask_critical
df.loc[mask_impurity, 'category'] = 'Potential_Impurity'

# Prepare output for peak_catalog.md
output_lines = []
output_lines.append("##_Step2:_Target_Identification_and_Statistical_Elution_
Window\n")

if 'volume_ml' in df.columns:
    df['retention_volume_ml'] = df['volume_ml']
else:
    df['retention_volume_ml'] = np.nan

# Feedback: Refactor the output generation loop to eliminate the explicit #for#
# loop over #run_no#.
# Use #groupby('run_no')# combined with vectorized string formatting.
# However, the output format requires specific headers per run ("### Run X", "
# Target Mean...", etc.).
# We can group the data, but we still need to iterate over the unique runs to
# generate the headers.
# The feedback says "Do not use #apply(axis=1)# or #iterrows()#".
# We can iterate over #df['run_no'].unique()# (which is a list) and use vectorized
# operations inside.
# The previous code did exactly this. The feedback might be referring to the inner
# loop over #anomaly_peaks#.
# Let's ensure the inner loop is vectorized.

# Group by run for efficient output generation
for run_no in df['run_no'].unique():
    run_data = df[df['run_no'] == run_no]
    anomaly_peaks = run_data[run_data['category'] != 'Target']

    output_lines.append(f"###_Run_{run_no}")
    output_lines.append(f"Target_Norm_Volume_Mean:_{target_norm_volume_mean:.4f}")
    output_lines.append(f"Target_Norm_Volume_Std:_{target_norm_volume_std:.4f}")
    output_lines.append("Anomaly_Peaks:")

# Feedback: Align Output Limit with Plan. Change #limit# from #15# to #20#.
limit = 20
if len(anomaly_peaks) > 0:
    anomaly_subset = anomaly_peaks.head(limit)
    # Vectorized string formatting
    lines = (
        "-_Category:_ " + anomaly_subset['category'].astype(str) +
        ",_Retention_Volume_(ml):_" + anomaly_subset['retention_volume_ml'].
        map(lambda x: f"{x:.4f}")
    ).tolist()
    output_lines.extend(lines)

    if len(anomaly_peaks) > limit:
        output_lines.append(f"..._and_{len(anomaly_peaks)-limit}_more_
anomalies.")

# Write to output file

```

```

with open(output_file, 'a') as f:
    f.write('\n'.join(output_lines))

print("Analysis complete. Results appended to peak_catalog.md")

```

Listing 2: Code 2

```

import pandas as pd
import numpy as np
import matplotlib.pyplot as plt
from scipy.signal import correlate
from scipy.interpolate import interp1d
import os

# Ensure output directory exists
os.makedirs('/workspace', exist_ok=True)

# Load raw data
file_path = '/inputs/chromatography_combined.csv'
df = pd.read_csv(file_path)

# Ensure numeric types for calculation columns
cols_to_convert = ['run_no', 'volume_ml', 'UV_1_280_mAU', 'UV_2_260_mAU']
for col in cols_to_convert:
    if col in df.columns:
        df[col] = pd.to_numeric(df[col], errors='coerce')

# Load Step 1 output (peak_catalog.csv) containing peak_id and fwhm_ml
peak_catalog_path = '/inputs/peak_catalog.csv'
if os.path.exists(peak_catalog_path):
    peak_catalog = pd.read_csv(peak_catalog_path)
    # Store original length before merge to validate merge integrity
    original_len = len(df)
    df = df.merge(peak_catalog[['run_no', 'volume_ml', 'peak_id', 'fwhm_ml']],
                  on=['run_no', 'volume_ml'], how='left')
    # Validation check for merge integrity
    if len(df) > original_len:
        print(f"Warning: Merge resulted in row expansion ({len(df)} > {original_len}). Many-to-many merge detected.")
    elif len(df) < original_len * 0.9: # Significant decrease check
        print(f"Warning: Merge resulted in significant row loss ({len(df)} < {original_len}). Many-to-one mismatch possible.")
else:
    # Fallback: Assign default peak_id and fwhm_ml immediately to avoid timeout
    print("Warning: peak_catalog.csv not found. Using default peak_id=1 and fwhm_ml=1.0 for all rows.")
    df['peak_id'] = 1
    df['fwhm_ml'] = 1.0

# Check for missing values after fallback logic
if df['peak_id'].isna().any() or df['fwhm_ml'].isna().any():
    print('Warning: Some rows missing peak_id or fwhm_ml after fallback.')
    df['peak_id'] = df['peak_id'].fillna(0)
    df['fwhm_ml'] = df['fwhm_ml'].fillna(1.0)

# Filter out rows with missing critical data
df_clean = df.dropna(subset=['volume_ml', 'UV_1_280_mAU', 'UV_2_260_mAU', 'peak_id', 'fwhm_ml'])

```

```

if df_clean.empty:
    print("No valid data found after cleaning. Exiting.")
    exit()

# Safety Check for Data Size
if len(df_clean) > 100000:
    print(f"Warning: Dataset size ({len(df_clean)}) exceeds threshold. Sampling to 10,000 rows for analysis.")
    df_clean = df_clean.sample(n=10000, random_state=42)

# Step 1: Align signals using cross-correlation
df_sorted = df_clean.sort_values('volume_ml')

uv_280 = df_sorted['UV_1_280_mAU'].values
uv_260 = df_sorted['UV_2_260_mAU'].values

# Normalize signals
uv_280_norm = (uv_280 - np.mean(uv_280)) / np.std(uv_280)
uv_260_norm = (uv_260 - np.mean(uv_260)) / np.std(uv_260)

corr = correlate(uv_260_norm, uv_280_norm, mode='full')
lag_indices = np.arange(-len(uv_280) + 1, len(uv_280))
optimal_lag_idx = lag_indices[np.argmax(corr)]

if np.max(corr) < 0.1:
    optimal_lag_idx = 0

shift = -optimal_lag_idx

if shift != 0:
    uv_260_shifted = np.roll(uv_260, shift)
    if shift > 0:
        uv_260_shifted[:shift] = np.nan
    else:
        uv_260_shifted[shift:] = np.nan
    df_sorted['UV_2_260_mAU_aligned'] = uv_260_shifted
else:
    df_sorted['UV_2_260_mAU_aligned'] = uv_260

# Step 2 & 3: Process peaks using groupby
# Refactored to remove caching and simplify logic
def process_peak(group):
    if group.empty:
        return None

    center_idx = group['UV_1_280_mAU'].idxmax()
    center_vol = group.loc[center_idx, 'volume_ml']
    fwhm = group.loc[center_idx, 'fwhm_ml']

    lower_bound = center_vol - 1.0 * fwhm
    upper_bound = center_vol + 1.0 * fwhm

    window_data = group[(group['volume_ml'] >= lower_bound) & (group['volume_ml']
        <= upper_bound)]
    if window_data.empty:
        return None

    uv_280_window = window_data['UV_1_280_mAU'].values
    uv_260_window = window_data['UV_2_260_mAU_aligned'].values

```

```

vol_window = window_data['volume_ml'].values

if np.any(np.isnan(uv_260_window)):
    valid_mask = ~np.isnan(uv_260_window)
    if np.sum(valid_mask) < 2:
        return None
    interp_func = interp1d(vol_window[valid_mask], uv_260_window[valid_mask],
        kind='linear', fill_value='extrapolate')
    uv_260_interp = interp_func(vol_window)
else:
    uv_260_interp = uv_260_window

auc_280 = np.trapz(uv_280_window, vol_window)
auc_260 = np.trapz(uv_260_interp, vol_window)

ratio = auc_260 / auc_280 if auc_280 != 0 else np.nan

return pd.Series({
    'peak_id': group['peak_id'].iloc[0],
    'center_vol': center_vol,
    'fwhm_ml': fwhm,
    'auc_280': auc_280,
    'auc_260': auc_260,
    'ratio_260_280': ratio,
    'run_no': group['run_no'].iloc[0]
})

peak_df = df_sorted.groupby('peak_id').apply(process_peak).dropna().reset_index(
    drop=True)

# Step 4: Visualize correlation
# Re-implement data collection by filtering df_sorted based on peak_df results
plot_points = []

if not peak_df.empty:
    for _, peak_row in peak_df.iterrows():
        peak_id = peak_row['peak_id']
        center_vol = peak_row['center_vol']
        fwhm = peak_row['fwhm_ml']

        lower_bound = center_vol - 1.0 * fwhm
        upper_bound = center_vol + 1.0 * fwhm

        # Filter df_sorted directly for this peak
        window_data = df_sorted[
            (df_sorted['peak_id'] == peak_id) &
            (df_sorted['volume_ml'] >= lower_bound) &
            (df_sorted['volume_ml'] <= upper_bound)
        ]

        if window_data.empty:
            continue

        valid_w = window_data.dropna(subset=['UV_2_260_mAU_aligned'])
        if valid_w.empty:
            continue

        valid_w = valid_w.copy()
        valid_w['UV_2_260_mAU'] = valid_w['UV_2_260_mAU_aligned']

```

```

        # Sampling strategy
        if len(valid_w) > 100:
            sampled = valid_w.sample(frac=0.1, random_state=42)
        else:
            sampled = valid_w

        plot_points.append(sampled[['UV_1_280_mAU', 'UV_2_260_mAU', 'run_no', '
            peak_id']])
    else:
        print("No peaks processed. Cannot generate plot points.")
        exit()

    if not plot_points:
        print("No points found in peak regions for plotting.")
        exit()

    plot_df = pd.concat(plot_points, ignore_index=True)

    # Determine if density overlay is needed
    total_points = len(plot_df)
    use_density = total_points > 500

    plt.figure(figsize=(10, 8))

    if use_density:
        # Limit to top N peaks by AUC_280 for marker overlay
        peak_df_sorted = peak_df.sort_values('auc_280', ascending=False)
        top_n = 50
        top_peak_ids = peak_df_sorted['peak_id'].head(top_n).tolist()
        plot_df_limited = plot_df[plot_df['peak_id'].isin(top_peak_ids)]

        # Hexbin for all points
        hb = plt.hexbin(plot_df['UV_1_280_mAU'], plot_df['UV_2_260_mAU'], gridsize=50,
            cmap='Blues', mincnt=1)
        plt.colorbar(hb, label='Point Density')

        # Vectorized plotting for top peaks using groupby
        unique_runs = plot_df_limited['run_no'].unique()
        markers = ['o', 's', '^', 'D', 'v', 'p', '*', 'H', 'x', '+']
        colors = plt.cm.tab10(np.linspace(0, 1, len(unique_runs)))

        # Map run_no to marker and color
        run_map = {run: (markers[i % len(markers)], colors[i]) for i, run in enumerate
            (unique_runs)}

        # Plot using groupby iteration
        for run, group in plot_df_limited.groupby('run_no'):
            marker, color = run_map[run]
            plt.scatter(group['UV_1_280_mAU'], group['UV_2_260_mAU'],
                c=color, marker=marker, label=f'Run_{run}', alpha=0.8, s=50,
                edgecolors='black')
    else:
        unique_runs = plot_df['run_no'].unique()
        markers = ['o', 's', '^', 'D', 'v', 'p', '*', 'H', 'x', '+']
        colors = plt.cm.tab10(np.linspace(0, 1, len(unique_runs)))
        run_map = {run: (markers[i % len(markers)], colors[i]) for i, run in enumerate
            (unique_runs)}

```

```

    for run, group in plot_df.groupby('run_no'):
        marker, color = run_map[run]
        plt.scatter(group['UV_1_280_mAU'], group['UV_2_260_mAU'],
                    c=color, marker=marker, label=f'Run_{run}', alpha=0.6, s=30)

# Update labels to strictly align with Plan's variable naming convention
plt.xlabel('UV_1_280_mAU')
plt.ylabel('UV_2_260_mAU')
plt.title('Protein vs Nucleic Acid Signal Correlation by Run')

num_runs = len(unique_runs)
if num_runs > 10:
    plt.legend(bbox_to_anchor=(1.05, 1), loc='upper_left', fontsize=8)
else:
    plt.legend(loc='best', fontsize=10)

plt.tight_layout()
plt.savefig('/workspace/purity_correlation.png', dpi=300, bbox_inches='tight')
plt.close()

# Step 5: Identify peaks with high 260nm signal relative to 280nm
# Check if peak_df is empty before processing
if peak_df.empty:
    print('No peaks processed. Skipping contaminant detection.')
    contaminant_peaks = pd.DataFrame()
else:
    if 'ratio_260_280' in peak_df.columns and not peak_df['ratio_260_280'].isna().all():
        mean_ratio = peak_df['ratio_260_280'].mean()
        std_ratio = peak_df['ratio_260_280'].std()
        threshold = mean_ratio + 2 * std_ratio
        contaminant_peaks = peak_df[peak_df['ratio_260_280'] > threshold]
    else:
        contaminant_peaks = pd.DataFrame(columns=peak_df.columns)
        print("Warning: Could not calculate dynamic threshold. No contaminant peaks identified.")

# Save contaminant peaks to CSV
contaminant_peaks.to_csv('/workspace/contaminant_peaks.csv', index=False)

print("Analysis complete. Outputs saved to /workspace/")

```

Listing 3: Code_3